







LYCOPENE OFFERS PROTECTION AGAINST OXIDATIVE DAMAGE IN FROZEN-THAWN BOVINE SEMEN

7th CASEE Conference
"The Role of Life Sciences
in Europe's 2020 Strategy"
May 22 - 24, 2016
Timisoara, Romania

Eva Tvrdá
Alica Mackovich
Hana Greifová
Norbert Lukáč

Department of Animal Physiology
Faculty of Biotechnology and Food
Sciences
Slovak University of Agriculture in Nitra
Slovakia

Sperm Cryopreservation

Benefits:


- Improvement of genetic quality and health condition of the animals
- Significant cost reduction for the production maintenance
- A valuable tool for assisted reproduction techniques in veterinary andrology
- An important technique for the protection of genetic potential and biodiversity endangered species
- Economic preservation of transgenic animal lines

ssues:

- Substantial loss of cellular mass (50-60%)
- Mechanical damage and thermal stress
- Oxidative stress



Oxidative Stress



Varicocele

Oxidative stress can damage this membrane preventing fertilisation.

Lycopene (LYC)

- The most abundant carotenoid found in nature, and a highly efficient antioxidant
- Proposed to be involved in processes related to the quenching of singlet oxygen and trapping of peroxy radicals
- Previous studies: LYC has the ability to improve sperm motility, membrane integrity and to protect DNA damage in spermatozoa





Other sources:














AIM OF THE STUDY

- o To evaluate the effects of lycopene on:
 - sperm motility
 - o reactive oxygen species (ROS) production
 - lipid peroxidation in cryopreserved bovine semen

Sample Collection and Processing

- 10 semen samples from 8 breeding bulls
- Minimum progressive motility 70 % and sperm concentration 1×10⁹ sperm/mL
- Two equal fractions diluted to 11x10⁶ sperm/mL
- Semen extender: Triladyl (egg yolk, TRIS, citric acid, sugar, buffers, glycerol and antibiotics) and diluted with distilled water
- Experimental group: the extender contained 1.5 mmol/L LYC
- Diluted semen samples \to 0.25 mL French straws \to 4 °C in 2 h \to freezing at a pre-programmed rate in a digital freezing machine










Spermatozoa Motility Analysis

Computer assisted semen analysis (CASA) system



Experimental procedure:

- 10 µL sample: Makler counting chamber
- 10 microscopic fields subjected to each analysis in order to include at least 300 cells
- Spermatozoa motility: motion > 5µm/s (%)

Reactive Oxygen Species (ROS) Generation

400 μL sample 10 μL DMSO

Luminol-based chemiluminiscence

15 repeated cycles of 1 min

Negative controls = sperm-free cryomedium/PBS + luminol

Positive controls = sperm-free cryomedium/PBS + H₂O₂ + luminol




Lipid Peroxidation Assessment




TBARS Method

- 100 µL sample + 4ml thiobarbituric acid (53%) dissolved in acetic acid (20%)
- Water bath: 1 hod. (100°C)
- Centrifugation: 8min/1800xg/4°C
- Plate reader: 530-540 nm


RESULTS: Spermatozoa Motility



RESULTS: ROS Production




RESULTS: Lipid Peroxidation



LYCOPENE EXHIBITS PROTECTIVE EFFECTS ON CRYOPRESERVED SPERMATOZOA

- Lower ROS generation
- Protection of membrane integrity
- Stimulation of spermatozoa motility



Potential practical application of the results

- Cryoconservation of mammalian, avian and reptilian spermatozoa
- Preservation of other somatic or reproductive cells and/or tissues
- Protection of genetic resources

THANK YOU FOR YOUR ATTENTION



This study was supported by the Scientific Grant Agency of the Ministry of Education of the Slovak Republic and of the Slovak Academy of Sciences VEGA Project no. 1/0857/14.